Operations for Sealing Film

SF-001-UC-100B

Product Description:

The Sealing Film is pressure sensitive covers, designed for use with 96-well and 384-well optical reaction plates only.

Feature:

Non-tacky to skin and gloves

High Clarity

Minimal to No Auto-fluorescence

DNase/RNase-free





Suitable for use between -70°C and 100°C (-94°F - 212°F)

Procedure:

- 1) Place the optical reaction plate in Plate Base, remove a single sealing film from the box.
- 2) Hold the sealing film back side up.
- 3) Fold back one of the solid cut line. Notice that the solid cut line is closer to the edge than the perforation.
- 4) Peel back the white protective backing from the sealing surface rapidly, taking care not to touch the sealing surface with your fingers.
- 5) With the sealing film positioned over the reaction plate slowly, holding the sealing film by the end-tabs.
- 6) Move the applicator slowly across the sealing film, while applying firm pressure to ensure good contact between the cover and the plate over the entire surface of the reaction plate.
- 7) While using the applicator to hold the edge of the cover in place, grasp the center of the end-tab and pull away sharply. Repeat for the other end-tab.
- 8) Transfer the sealed plate to the instrument.
- 9) Run the thermal cycling protocol.

Note:

- 1) Always handle the optical adhesive cover by the end-tabs.
- 2) Improper peeling may result in haziness, but will not affect results. Haziness will disappear when it comes into contact with the heated cover in the instrument.
- 3) This product is pressure sensitive covers, which require the application of pressure by the user to ensure a tight, evaporation-free seal.

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